

Congratulations!

You have just received the new Zoom Chairside Whitening Procedure kit.

Discus is proud to present the new and improved Zoom Chairside Whitening Procedure kit. This redesigned kit makes the Zoom procedure easier to administer and safer for your patients.

Based on the input of hundreds of Zoom doctors like yourself, we are confident you will find this new kit an even better value than before.



For training on the new kit, visit zoomtraining.com

ZOOM!

DISCUS DENTAL®

© 2009 Discus Dental, LLC. All rights reserved.
To be dispensed by or on the order of a dental professional only

We've made things easier, faster and safer by:

- 1 Including a stylish new Patient Post Care & Maintenance kit, making it easy for patients to manage sensitivity and maintain their white smile.
- 2 Engineering a new IsoPrep™ retractor designed to provide complete lip coverage and easier access.
- 3 Reformulating the whitening gel to optimize pH, eliminating the need for pre-treatment swabs.
- 4 Including simple step-by-step pictorial instructions on isolation and whitening.
- 5 Adding a surgical suction tip for easier removal of the Zoom whitening gel without disturbing Liquidam.
- 6 Adding cotton rolls so there's no need to roll gauze for placement in the upper and lower vestibules.
- 7 Reorganizing the kit contents in procedural order.
- 8 Introducing a new Liquidam with improved viscosity and 30% – 40% more UV protection.

